
Recombinant human IgG3 Fc protein,(HEK293)

Catalog Number: bs-43577P

Concentration: >1mg/ml

Species: Human

AA Seq: 99-377/377

Predicted MW: 31.2

Tags: Tag Free

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH=7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.