

---

## Recombinant human IgG4 Fc protein,(HEK293)

Catalog Number: bs-43578P

Concentration: >1mg/ml

AA Seq: 99-327/327

Predicted MW: 25.8

Tags: Tag Free

Purity: >90% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: PBS (pH=7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.