

---

## Recombinant ASFV p30 protein, N-His

Catalog Number: bs-41382P

Concentration: >0.5 mg/ml

AA Seq: 1-204/204

Predicted MW: 23.5

Detected MW: 33 kDa

Tags: N-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** Modifies the subcellular distribution of heterogeneous nuclear ribonucleoprotein K (HNRNPK) and may contribute to modulate HNRNPK functions related to processing and export of mRNAs during ASFV infection.

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.