
Recombinant SARS-Cov-2 Spike RBD protein, His, FITC conjugated

Catalog Number: bs-41516P-FITC

AA Seq: 319-537/1273

Predicted MW: 24.7

Tags: His tag

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 0.2µm filtered solution of 0.01M PBS (pH7.4) with 0.02% Proclin300.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: The SARS-CoV-2 spike (S) protein is the target of vaccine design efforts to end the COVID-19 pandemic. Despite a low mutation rate, isolates with the D614G substitution in the S protein appeared early during the pandemic, and are now the dominant form worldwide. Here, we analyze the D614G mutation in the context of a soluble S ectodomain construct.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.