

---

## Recombinant human Serpin B4 protein, His (HEK293)

Catalog Number: bs-43150P

Concentration: >1mg/ml

AA Seq: 2-390/390

Predicted MW: 46.8

Tags: C-His

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: May act as a protease inhibitor to modulate the host immune response against tumor cells.

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.