
Recombinant human CD28 protein, C-hFc (HEK193)

Catalog Number: bs-43587P

Concentration: >1mg/ml

Species: Human

AA Seq: 19-152/220

Predicted MW: 41.4

Tags: C-hFc

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2011]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.