

bs-4581R

[Primary Antibody]

Listeria tropina Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit Clonality: Polyclonal Target: Listeria tropina Purification: affinity purified by Protein A Concentration: 1mg/ml Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	Isotype: IgG	Applications: ELISA (1:5000-10000) Reactivity: (predicted: Listeria)
---	---------------------	---

— SELECTED CITATIONS —

- **[IF=2.6]** Hong Hailong. et al. Listeria monocytogenes: possible mechanism of infection of goat uterus and its effects on uterine autophagy and cell apoptosis. FRONT VET SCI. 2024 Aug;11: IHC ;Goat. 39220769