

---

## Human IL-6 Mouse mAb

Cat. No:	V0401
Name:	Human IL-6 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	6F5K
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA
Antigen:	Recombinant human IL-6 (C-6x His-Tag) ( <a href="#">bs-41177P</a> )