

---

## AT-2 Peptide

Catalog Number: bs-0587P

AA Seq: (DRVYIHPF-GG)8K4K2KG

Activity: Yes

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: bs-0587P is a Eight branched multiple antigenic peptide of Angiotensin II.

### PRODUCT SPECIFIC PUBLICATIONS

---

**[IF=1.96]** Dong, Fang, et al. "Amperometric Immunosensor Based on Carbon Nanotubes/Chitosan Film Modified Electrodes for Detection of Human Leptin." International Journal of Electrochemical Science 9 (2014): 6924-6935. Other ; ="" . 10.1016/S1452-3981(23)10941-2