

**bs-19446R**

**[ Primary Antibody ]**

## PRTFDC1 Rabbit pAb

**Bioss**  
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn


techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> WB (1:500-2000)  <b>Reactivity:</b> Rat (predicted: Human, Rabbit, Pig, Dog)  <b>Predicted MW.:</b> 26 kDa  <b>Subcellular Location:</b> Cytoplasm
<b>Clonality:</b> Polyclonal		
<b>GeneID:</b> 56952	<b>SWISS:</b> Q9NRG1	
<b>Target:</b> PRTFDC1		
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from human PRTFDC1: 151-225/225.		
<b>Purification:</b> affinity purified by Protein A		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

### — VALIDATION IMAGES —



Sample: Testis(Rat) Cell Lysate at 40 ug Primary:  
Anti-PRTFDC1 (bs-19446R) at 1/300 dilution  
Secondary: IRDye800CW Goat Anti-Rabbit IgG at  
1/20000 dilution Predicted band size: 26 kD  
Observed band size: 21 kD