

bs-19447R

[Primary Antibody]

PRTG Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 283659</p> <p>Target: PRTG</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human PRTG: 311-410/1150. < Extracellular ></p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: This gene encodes a member of the immunoglobulin superfamily. The encoded transmembrane protein has been associated with the development of various tissues, especially neurogenesis. It has been suggested that this gene may be associated with attention deficit hyperactivity disorder (ADHD). [provided by RefSeq, Nov 2014]</p>	<p>Isotype: IgG</p> <p>SWISS: Q2VWP7</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Sheep, Cow, Chicken, Dog, Horse)</p> <p>Predicted MW.: 124 kDa</p> <p>Subcellular Location: Cell membrane ,Cytoplasm</p>
---	--	---