
PRTG Rabbit pAb

Catalog Number: bs-19447R

Target Protein: PRTG

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse, Rat, Sheep, Cow, Chicken, Dog, Horse)

Predicted MW: 124 kDa

Entrez Gene: 283659

Swiss Prot: Q2VWP7

Source: KLH conjugated synthetic peptide derived from human PRTG: 311-410/1150.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a member of the immunoglobulin superfamily. The encoded transmembrane protein has been associated with the development of various tissues, especially neurogenesis. It has been suggested that this gene may be associated with attention deficit hyperactivity disorder (ADHD). [provided by RefSeq, Nov 2014]