

Recombinant SARS-Cov-2 N protein (del204, del215), His

Catalog Number: bs-41491P

Concentration: >0.5 mg/ml

AA Seq: 1-419/419

Predicted MW: 46

Detected MW: 50 kDa

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC


Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: The SARS-CoV-2 spike (S) protein is the target of vaccine design efforts to end the COVID-19 pandemic. Despite a low mutation rate, isolates with the D614G substitution in the S protein appeared early during the pandemic, and are now the dominant form worldwide. Here, we analyze the D614G mutation in the context of a soluble S ectodomain construct.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.