
Recombinant Outer Membrane protein, Trx & S & His

Catalog Number: bs-49105P

Concentration: >0.5 mg/ml

AA Seq: 1-231/231

Predicted MW: 39.8

Tags: Trx & S & His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 20mM Tris-HCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.