

---

## **BOLA2 Antibody Blocking Peptide**

Catalog Number: bs-8422P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene is located within a region of a segmental duplication on chromosome 16p. The product of this gene belongs to the eukaryotic subfamily of the Bola-like proteins. This gene encodes the Bola-like protein 2. The Bola-like proteins are widely conserved from prokaryotes to eukaryotes, and these proteins seem to be involved in cell proliferation or cell-cycle regulation. This duplicated region contains a complete copy of the short splice variant of the BOLA2 gene and a partial copy of the long splice variant of the BOLA2 gene.