
UNC93B Antibody Blocking Peptide

Catalog Number: bs-8504P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: UNC93B is a protein with similarity to the *C. elegans* unc93 protein. The Unc93 protein is involved in the regulation or coordination of muscle contraction in the worm. Defects in UNC93B are associated with susceptibility to herpes simplex encephalitis (HSE) which is a rare complication of human herpesvirus 1 (HHV-1) infection.