

---

## Vapreotide, Acetate Rabbit pAb

Catalog Number: bs-37050R

Target Protein: Vapreotide, Acetate

Concentration: 1mg/1ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Vapreotide)

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Vapreotide (Sanvar) is a synthetic somatostatin analog. It is used in the treatment of esophageal variceal bleeding in patients with cirrhotic liver disease and AIDS-related diarrhea.