

**bs-37050R**

**[ Primary Antibody ]**

## Vapreotide, Acetate Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Vapreotide, Acetate

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/1ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** Vapreotide (Sanvar) is a synthetic somatostatin analog. It is used in the treatment of esophageal variceal bleeding in patients with cirrhotic liver disease and AIDS-related diarrhea.

**Applications:** IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

**Reactivity:** (predicted: Vapreotide)