
VOPP1 Rabbit pAb

Catalog Number: bs-12760R

Target Protein: VOPP1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Sheep, Cow, Dog, Horse)

Predicted MW: 17 kDa

Entrez Gene: 81552

Swiss Prot: Q96AW1

Source: KLH conjugated synthetic peptide derived from human VOPP1: 41-140/172.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: VOPP1 (Vesicular, Overexpressed In Cancer, Prosurvival Protein 1) is a Protein Coding gene.

GO annotations related to this gene include signal transducer activity.