
Normal Guinea pig Serum

Catalog Number: C-0358

Form: Liquid

Applications: Isotype Control, Blocking Assay, etc.

Storage: Neat serum with 0.02%Proclin300.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This serum is a pooled sample collected from healthy adult animals. The serum is heat-treated and centrifuged to remove precipitates and then filtered. Each serum is tested with the appropriate antibody to check for possible cross-reactivities. The sera are supplied undiluted with 0.08% sodium azide as a preservative. When diluted, this reagent can be used for blocking non-specific binding or as an antibody diluent.