

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Normal Chicken Serum

Catalog Number: C-0310

Form: Lyophilized

Applications: Isotype Control, Blocking Assay, etc.

Storage: Neat serum with 0.02%Proclin300.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Normal chicken serum. Whole serum is obtained from non-immunized (preimmune) healthy

chicken. The serum is clarified through a 0.45 micron low protein binding filter.