
Human IL-31 ELISA Kit

产品编号: bsk11076

种属: Human

线性范围: 125 - 8,000 pg/mL

应用范围: S/P/CC

检测限: 60 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).

产品介绍: IL31, which is made principally by activated Th2-type T cells, interacts with a heterodimeric receptor consisting of IL31RA (MIM 609510) and OSMR (MIM 601743) that is constitutively expressed on epithelial cells and keratinocytes. IL31 may be involved in the promotion of allergic skin disorders and in regulating other allergic diseases, such as asthma (Dillon et al., 2004 [PubMed 15184896]).[supplied by OMIM, Mar 2008]