
Human IL-34 ELISA Kit

产品编号: bsk11078

种属: Human

线性范围: 64 - 4,000 pg/mL

应用范围: S/P/CC

检测限: 30 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles
(Shipped with wet ice.).

产品介绍: Interleukin-34 is a cytokine that promotes the differentiation and viability of monocytes and macrophages through the colony-stimulating factor-1 receptor (CSF1R; MIM 164770) (Lin et al., 2008 [PubMed 18467591]).[supplied by OMIM, May 2008]