
Human PDGF-BB ELISA Kit

产品编号: bsk11090

种属: Human

线性范围: 31.20 - 2,000 pg/mL

应用范围: S/P/CC

检测限: 15 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).

产品介绍: Platelet Derived Growth Factor (PDGF) is the major serum mitogen for cells of mesenchymal origin in humans. The biologically active protein is a dimer composed of two related polypeptides designated A and B, the dimers being connected by disulfide bonds. The PDGF protein has been implicated both directly as well as indirectly in several pathological states including neoplasia, arthritis, arteriosclerosis and bone marrow sclerosis.

PRODUCT SPECIFIC PUBLICATIONS

[IF=7.31] Panshi Jin. et al. Platelets Facilitate Wound Healing by Mitochondrial Transfer and Reducing Oxidative Stress in Endothelial Cells. OXID MED CELL LONGEV. 2023;2023:2345279 ELISA ; Human . 36860732