
Rat ICAM-1/CD54 ELISA Kit

产品编号: bsk13022

种属: Rat

线性范围: 62.50 - 4,000 pg/mL

应用范围: S/P/CC

检测限: 31 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles
(Shipped with wet ice.).

产品介绍: This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor. [provided by RefSeq, Jul 2008].