

---

## Human CCL1/I-309 ELISA Kit

产品编号: bsk11150

种属: Human

线性范围: 15.60 - 1,000 pg/mL

应用范围: S/P/CC

检测限: 7 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).

产品介绍: This gene is one of several Cys-Cys (CC) cytokine genes clustered on the q-arm of chromosome 17. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The cytokine encoded by this gene displays chemotactic activity for eosinophils, but not mononuclear cells or neutrophils. This eosinophil specific chemokine assumed to be involved in eosinophilic inflammatory diseases such as atopic dermatitis, allergic rhinitis, asthma and parasitic infections. [provided by RefSeq, Jul 2008].