
Human IL-18 BPα ELISA Kit

产品编号: bsk11171

种属: Human

线性范围: 93.75-6000 pg/mL

应用范围: S/P/CC

检测限: 47 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).

产品介绍: IL18 binding protein binds to IL18, prevents the binding of IL18 to its receptor, and thus inhibits IL18 induced IFN gamma production. It is constitutively expressed and secreted in mononuclear cells and can be enhanced by IFN gamma. An elevated level of this protein is detected in the intestinal tissues of patients with Crohn disease. Four alternatively spliced transcript variants encoding distinct isoforms have been reported.