

---

## Mouse BAFF/BLyS/TNFSF13B ELISA Kit

产品编号: bsk12049

种属: Mouse

线性范围: 62.50-4000 pg/mL

应用范围: S/P/CC

检测限: 30 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles  
(Shipped with wet ice.).

产品介绍: The protein encoded by this gene is a cytokine that belongs to the tumor necrosis factor (TNF) ligand family. This cytokine is a ligand for receptors TNFRSF13B/TACI, TNFRSF17/BCMA, and TNFRSF13C/BAFFR. This cytokine is expressed in B cell lineage cells, and acts as a potent B cell activator. It has been also shown to play an important role in the proliferation and differentiation of B cells. [provided by RefSeq].

### PRODUCT SPECIFIC PUBLICATIONS

---

[IF=5.195] Xiaofeng Niu. et al. Panlongqi tablet suppresses adjuvant-induced rheumatoid arthritis by inhibiting the inflammatory response in vivo and in vitro. J ETHNOPHARMACOL. 2023 Feb;;116250 ELISA ; Mouse . 36791928