

bs-0493R

[Primary Antibody]

SPATA3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 130560</p> <p>Target: SPATA3</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SPATA3: 101-183/183.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: TSARG3 may be involved in inhibiting testis spermatogenesis apoptosis.</p>	<p>Isotype: IgG</p> <p>SWISS: Q8NHX4</p>	<p>Applications: WB (1:500-2000) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit)</p> <p>Predicted MW.: 20 kDa</p> <p>Subcellular Location: Nucleus</p>
---	--	---