## bs-0493R

## [ Primary Antibody ]

## BIOSS ANTIBODIES

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## **SPATA3 Rabbit pAb**

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GeneID: 130560 SWISS: Q8NHX4

Target: SPATA3

**Immunogen:** KLH conjugated synthetic peptide derived from human SPATA3:

101-183/183.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

**Background:** TSARG3 may be involved in inhibiting testis spermatogenesis

apoptosis.

**Applications: WB** (1:500-2000)

**ELISA** (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Rabbit)

Predicted MW.: 20 kDa

Subcellular Nucleus