

---

## Mouse NGAL/Lipocalin-2 ELISA Kit

产品编号: bsk12054

种属: Mouse

线性范围: 7.8-500 pg/mL

应用范围: S/P/CC

检测限: 4 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).

产品介绍: Lipocalin 2 is a member of the lipocalin family which encompass more than 25 members (including proteins like retinol binding protein,  $\alpha 1$  microglobulin,  $\beta$  lactoglobulin, apolipoprotein D and odorant binding protein). The members of the lipocalin family are characterized by their ability to bind small lipophilic substances in their hydrophobic core, and thereby serve as transporters of substances.