
Mouse IFN-beta ELISA Kit

产品编号: bsk12073

种属: Mouse

线性范围: 12.5 - 800 pg/mL

应用范围: S/P/CC

检测限: 6 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles
(Shipped with wet ice.).

产品介绍: Interferon-Beta, IFN-Beta mainly by the fibroblast cell factor, belonging to a single gene encodes a protein that can interfere with the virus RNA or DNA replication and inhibit virus growth, and significantly enhance NK cell activity, through the promotion of MHC I molecule expression and enhance the CTL of the virus-infected cells role in the identification and destruction.

PRODUCT SPECIFIC PUBLICATIONS

[IF=14.3] Yang Zhou. et al. T Cell-Derived Apoptotic Extracellular Vesicles Hydrolyze cGAMP to Alleviate Radiation Enteritis via Surface Enzyme ENPP1. ADV SCI. 2024 Jun;;2401634 ELISA ; Mouse . 38888507

[IF=6.8] Hong-Fei Wang. et al. Jujuboside B alleviates acetaminophen-induced hepatotoxicity in mice by regulating Nrf2-STING signaling pathway. ECOTOX ENVIRON SAFE. 2024 Jan;269:115810 ELISA ; Mouse . 38100849