

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Mouse IFN-beta ELISA Kit

产品编号: bsk12073

种属: Mouse

线性范围: 12.5 - 800 pg/mL

应用范围: S/P/CC

检测限: 6 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles

(Shipped with wet ice.).

产品介绍: Interferon-Beta, IFN-Beta mainly by the fibroblast cell factor, belonging to a single gene

encodes a protein that can interfere with the virus RNA or DNA replication and inhibit virus

growth, and significantly enhance NK cell activity, through the promotion of MHC I

molecule expression and enhance the CTL of the virus-infected cells role in the identification

and destruction.

PRODUCT SPECIFIC PUBLICATIONS

[IF=14.3] Yang Zhou. et al. T Cell-Derived Apoptotic Extracellular Vesicles Hydrolyze cGAMP to Alleviate Radiation Enteritis via Surface Enzyme ENPP1. ADV SCI. 2024 Jun;:2401634 ELISA; Mouse . 38888507

[IF=14.1] Zhengyan Wang. et al. Cellular Adaptation to Mechanical Stress Emerges via Cell Shrinkage Triggered by Nonlinear Calcium Elevation. ADV SCI. 2025 Jul;:e03659 ELISA; Mouse . 40679358

[IF=6.8] Hong-Fei Wang. et al. Jujuboside B alleviates acetaminophen-induced hepatotoxicity in mice by regulating Nrf2-STING signaling pathway. ECOTOX ENVIRON SAFE. 2024 Jan;269:115810 ELISA; Mouse . 38100849