

---

## Mouse Thrombopoietin ELISA Kit

产品编号: bsk12082

种属: Mouse

线性范围: 62.5 - 4000 pg/mL

应用范围: S/P/CC

检测限: 32 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).

产品介绍: Megakaryocytopoiesis is the cellular development process that leads to platelet production. The protein encoded by this gene is a humoral growth factor that is necessary for megakaryocyte proliferation and maturation, as well as for thrombopoiesis. This protein is the ligand for MLP/C\_MPL, the product of myeloproliferative leukemia virus oncogene. Mutations in this gene are the cause of thrombocythemia 1. Alternate splicing results in multiple transcript variants of this gene.[provided by RefSeq, Jun 2012].