
Mouse CXCL1 ELISA Kit

产品编号: BSKM1005

种属: Mouse

线性范围: 15.6-1000 pg/mL

应用范围: S/P/CC

检测限: 7 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles
(Shipped with wet ice.).

PRODUCT SPECIFIC PUBLICATIONS

[IF=6.208] Xuena Yang. et al. A CpG-Oligodeoxynucleotide Suppresses Th2/Th17 Inflammation by Inhibiting IL-33/ST2 Signaling in Mice from a Model of Adoptive Dendritic Cell Transfer of Smoke-Induced Asthma. INT J MOL SCI. 2023 Jan;24(4):3130 ELISA ; Mouse . 36834541

[IF=5.1] Liu Jing. et al. Interplay of IL-33 and IL-35 Modulates Th2/Th17 Responses in Cigarette Smoke Exposure HDM-Induced Asthma. INFLAMMATION. 2023 Sep;:1-18 ELISA ; Mouse . 37737467