
Human T3 ELISA Kit

产品编号: bsk31006

种属: Human

线性范围: 0—2000ng/dl

应用范围: 人血清

检测限: 25ng/dl

适用样品基质: 人血清

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Shipped with wet ice.).

产品介绍: T3(Triiodothyroxine)although only a minority of the total thyroxine, but the T3 is very important to check thyroid function. In the biochemical function of T4 than strong, and in the case of hyperthyroidism, T3 concentration in the blood are often marked increase (or even up to 5 times the normal).