

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## **Human IL-22 ELISA Kit**

产品编号: bsk11030

种属: Human

线性范围: 15.6 - 1,000 pg/mL

应用范围: S/P/CC

检测限: 6.0 pg/mL

适用样品基质: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological

fluids.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles

(Shipped with wet ice.).

产品介绍: This gene was identified by its overexpression specifically in herpesvirus samimiri-

transformed T cells. The encoded protein is a member of the IL10 family of cytokines. It is a secreted protein and may function as a homodimer. This protein is thought to contribute to the transformed phenotype of T cells after infection by herpesvirus samimiri. [provided by

RefSeq].