

---

## Human IL-24/MDA-7 ELISA Kit

产品编号: bsk11060

种属: Human

线性范围: 78.1-5000pg/mL

检测限: 27pg/mL

适用样品基质: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

保存条件: Store at 4°C & -20°C for 12 months.

产品介绍: This gene encodes a member of the IL10 family of cytokines. It was identified as a gene induced during terminal differentiation in melanoma cells. The protein encoded by this gene can induce apoptosis selectively in various cancer cells. Overexpression of this gene leads to elevated expression of several GADD family genes, which correlates with the induction of apoptosis. The phosphorylation of mitogen-activated protein kinase 14 (MAPK7/P38), and heat shock 27kDa protein 1 (HSPB2/HSP27) are found to be induced by this gene in melanoma cells, but not in normal immortal melanocytes. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]