
Human IL-17C ELISA Kit

产品编号: bsk11054

种属: Human

线性范围: 125.00 - 8,000 pg/mL

应用范围: S/P/CC

检测限: 60 pg/mL

适用样品基质: cell culture supernates, serum, and plasma.

保存条件: Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles
(Shipped with wet ice.).

产品介绍: The protein encoded by this gene is a T cell-derived cytokine that shares the sequence similarity with IL17. This cytokine was reported to stimulate the release of tumor necrosis factor alpha and interleukin 1 beta from a monocytic cell line. The expression of this cytokine was found to be restricted to activated T cells. [provided by RefSeq].