
MRO Rabbit pAb

Catalog Number: bs-17764R

Target Protein: MRO

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human)

Predicted MW: 29 kDa

Entrez Gene: 83876

Swiss Prot: Q9BYG7

Source: KLH conjugated synthetic peptide derived from human MRO: 1-100/248.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene is specifically transcribed in males before and after differentiation of testis, and the encoded protein may play an important role in a mammalian sex determination.

Multiple transcript variants encoding different isoforms have been found for this gene.

[provided by RefSeq, Jul 2008]