

bs-7042R

[Primary Antibody]

ELF1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 1997</p> <p>Target: ELF1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human ELF1: 201-300/619.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: E74 like Factor 1 is a transcription factor that may be required for the T-cell receptor mediated activation of HIV2 gene expression. It binds to two purine-rich motifs in the HIV-2 enhancer and transcriptionally activates the LYN and mouse BLK promoters.</p>	<p>Isotype: IgG</p> <p>SWISS: P32519</p> <p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: Rat (predicted: Human, Mouse, Pig, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 67 kDa</p> <p>Subcellular Location: Nucleus</p>
--	--