

---

## CMG1 Rabbit pAb

Catalog Number: bs-7050R

Target Protein: CMG1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Cow, Horse)

Predicted MW: 69 kDa

Entrez Gene: 80173

Source: KLH conjugated synthetic peptide derived from human CMG1: 51-150/600.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: No data available.