
Malic Dehydrogenase from porcine heart

产品编号: D10113

CAS: 9001-64-3

保存条件: Store at 2-8°C.

产品介绍: CAS号: 9001-64-3

分子量: 34KDa (SDS-PAGE)

MDL: MFCD00011351

抑制剂: Ag⁺, Hg²⁺

外观/性状: Suspension in 2.8 M (NH₄)₂SO₄ solution, pH 6.0

activity: ≥600 units/mg protein (biuret), ammonium sulfate suspension Application:

Malic Dehydrogenase from porcine heart has been used:

- in qualitative protein binding measurements
- to test internally calibrated electrochemical continuous enzyme assay (ICECEA) with model enzyme pair
- to investigate the effect of chaperone on the refolding of heat-denatured malate dehydrogenase