
Biotin-NHS

产品编号: D-9811

CAS: 35013-72-0

保存条件: Store at -20°C, protect from light.

产品介绍: May be used in the biotinylation of proteins and peptides. Typically coupled to primary amine in the pH range 6.5-8.5.

PRODUCT SPECIFIC PUBLICATIONS

[IF=8.758] Lei Bi. et al. Microwave-Assisted Synthesis of Hollow Microspheres with Multicomponent Nanocores for Heavy-Metal Removal and Magnetic Sensing. Acs Appl Mater Inter. 2020;12(41):46779–46787
Other ; . 32965095