
Sorafenib

产品编号: D10317

CAS: 284461-73-0

保存条件: Store at 2-8°C.

产品介绍: 纯度: ≥98%

CAS号: 284461-73-0

分子式: C₂₁H₁₆ClF₃N₄O₃

分子量: 464.83

外观/性状: 白色至灰白色固体

敏感性: 易吸潮

溶解性: Soluble in DMSO (~20 mg/ml), Hot Methanol, DMF (~20 mg/ml), 1:2 solution of DMSO:PBS(PH 7.2) (~0.3 mg/ml), 100% ethanol.

用途: An inhibitor of tyrosine protein kinases that targets the Raf/Mek/Erk pathway.