
Cardiac Troponin I (cTnI) Mouse mAb, HRP conjugated

Cat. No:	V3201-HRP
Name:	Cardiac Troponin I (cTnI) Mouse mAb, HRP conjugated
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	17C8L
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Concentration:	1mg/ml
Buffer:	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Antigen:	Recombinant human Cardiac Troponin I (cTnI)