

---

## HT培养基添加剂 (50×) HT Media Supplement (50×)

产品编号: C7092

保存条件: Store at -20°C.

产品介绍: Application

HT Media supplement is a cocktail hypoxanthine and thymidine qualified for use in hybridoma selection and cell culture applications. HT Media Supplement (50×) Hybri-Max™ has been used as a supplement in INS-1 β-cells co-culture medium. It is suitable for use as medium for hybridoma cells and as Dulbecco's Modified Eagle Medium (DMEM) for HeLa S3 and MTH1-KO cells.

### Reconstitution

Reconstitute contents of vial with 10 ml sterile cell culture medium. Stock solution is sufficient to prepare 500 ml medium. Final working concentration: 100 μM hypoxanthine, 16 μM thymidine.