

bsm-63121R

[Primary Antibody]

alpha/beta Synuclein Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 1C4

GeneID: 6622

SWISS: P37840

Target: alpha/beta Synuclein

Immunogen: A synthesized peptide derived from human alpha Synuclein: 1-28.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively.

Applications: **WB** (1:1000-1:2000)

Flow-Cyt (1:20-1:100)

ICC/IF (1:50-1:200)

Reactivity: Human, Mouse, Rat

Predicted
MW.: 14 kDa

Subcellular Cell membrane ,Cytoplasm
Location: ,Nucleus