

bsm-63138R

[Primary Antibody]

MMP8 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

CloneNo.: 10H7

GeneID: 4317

SWISS: P22894

Target: MMP8

Immunogen: A synthesized peptide derived from human MMP8: 150-200.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

Background: MMP8 Can degrade fibrillar type I, II, and III collagens. Belongs to the peptidase M10A family.

Applications: **WB** (1:1000-1:2000)

IHC-P (1:100-1:200)

IHC-F (1:100-1:200)

IF (1:50-1:200)

Flow-Cyt (1:20-1:100)

ICC/IF (1:50-1:200)

IP (1:20-1:50)

Reactivity: Human, Mouse, Rat

Predicted MW.: 53 kDa

Subcellular Location: Secreted ,Extracellular matrix ,Cytoplasm