

**bsm-63185R**

**[ Primary Antibody ]**

## RAF1 Recombinant Rabbit mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Recombinant

**CloneNo.:** 6B6

**GeneID:** 5894

**SWISS:** P04049

**Target:** RAF1

**Immunogen:** A synthesized peptide derived from human c Raf: 1-67.

**Purification:** affinity purified by Protein A

**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.

**Background:** Serine/threonine-protein kinase that acts as a regulatory link between the membrane-associated Ras GTPases and the MAPK/ERK cascade, and this critical regulatory link functions as a switch determining cell fate decisions including proliferation, differentiation, apoptosis, survival and oncogenic transformation.

**Applications:** **WB** (1:500-1:2000)  
**ICC/IF** (1:50-1:200)

**Reactivity:** Human

**Predicted**  
**MW.:** 73 kDa

**Subcellular** Cell membrane ,Cytoplasm  
**Location:** ,Nucleus